





Molar Concentration

Consider the following:



1.0 L of H₂O



0.50 L of H₂O

a) How many moles are in 1.0 L of H_2O ?

Answer: 1 mole per liter

b) Now, how many moles are in 1.0 L of H_2O ?

Answer: 2 moles per liter

c) How many moles are in 0.5 L of H_2O ?

Answer: 2 moles per 0.5 L

Or 4 moles per liter



Molar Concentration

Mole formula #3

$$C = \frac{n}{V}$$
 or $n = C \times V$
 $n = \# \text{ of moles (mol)}$
 $v = \text{ volume (L)}$

where

 $C = \text{concentration (mol/L or } \mathbf{M})$

- [] square bracket also a symbol of concentration
 - e.g. [HCl_(aa)] means "the concentration of HCl....."
 - See sample problem #4 and #5 on page 283-284



Sample problem

- A sample of ammonia (NH₃) solution has a concentration of 2.6 mol/L.
 - a) How many moles of ammonia would be in a 250 mL bottle?
 - b) What mass of ammonia was used?

```
Given: C = 2.6 \text{ mol/L}, V = 250 \text{ mL} = 0.250 \text{ L}

n = C \times V

= 2.6 \text{ mol/L} \times 0.250 \text{ L}

= 0.65 \text{ mol}

m (NH_3) = n \times M

= 0.65 \text{ mol} \times 17.0 \text{ g/mol}

= 11.0 \text{ g}

= 1.1 \times 10^1 \text{ g} (2 \text{ sig digits})
```



Your Turn!

- Read through 281-290 as review of lesson
 - See Summary chart on pg 290
- Q#1-6, 8 on pg 284
- Q#12,13, 15,16 on pg 287
- Q#19, 20 on pg 290



Warm-up Question

□ Try Q#10 on page 291



Preparing Solutions – from a Solid

Creating a Solution from a Solid

- How would you create 500 mL of a 2.5 mol/L solution of sodium hydroxide?
- 1. Calculate the number of moles in the solution
 - $n = C \times V$
 - $= 2.5 \text{ mol/L} \times 0.500 \text{ L}$
 - = 1.25 mol
- 2. The number of moles can not be measured?
 - i. What can we measure? I mass
 - ii. Find the mass of solid required
 - $m = n \times M$
 - $= 1.25 \text{ mol} \times 40.0 \text{ g/mol}$
 - = 50 g of NaOH (2 sig digits)

Use a balance to mass out this amount of solid NaOH.



3. Obtain a 500 mL volumetric flask.

Add the measured amount of solid to the flask

50.0 g




Add approx 250 mL of distilled water to the flask. Cap and invert to mix (several times)

5. Add remaining distilled water to the <u>fill line</u>. Cap and invert again.



Dilutions of a Stock solution to form a weaker solution




10 mL of 0.25 M HCl

$$n = C \times V$$

 $= 0.25 \text{ mol/L} \times 0.010 \text{ L}$

$$= 2.5 \times 10^{-3} \text{ mol}$$



DILUTE by adding 90 mL of water for a TOTAL volume of 100 mL

$$C = n / V$$

= 0.0025 mol / 0.100 L

= 0.025 mol/L

 $= 2.5 \times 10^{-2} \text{ M}$

Note: After adding more water, we still have the same number of moles, but conc. decreases due to increase in volume



moles before = moles after

$$egin{aligned} egin{aligned} egin{aligned\\ egin{aligned} egi$$

(This formula can only be used for same solution on both sides)

- Sample Problem: You have a solution of 6.0 M HCl and you want to create 1.0 L of a 1.5 M acid solution for Mr. Gladden.
- Given: $C_1 = 6.0 \text{ mol/L}$ $C_2 = 1.5 \text{ mol/L}$ $V_1 = ? L$ $V_2 = 1.0 \text{ L}$

Equation

$$C_1V_1 = C_2V_2$$

$$V_1 = \frac{(1.5 \text{ mol/L})(1.0 \text{L})}{6.0 \text{ mol/L}}$$

- = 0.250 L
- $= 250 \, mL$
- you need 250 mL of 6.0 M HCl and 750 mL of water.



Method

- Obtain a 1000 mL volumetric flask. Add about 500 mL of distilled water.
- Measure 250 mL of concentrated HCl_(aq) with a graduated solution. Add carefully to the flask.
- Stopper, swirl and invert to mix solute and solvent
- Add more distilled water to fill line. Stopper, swirl and mix.

See sample problem #2 on pg 304.

WFT: Q#1,2,4,5a on pg 302

- Q#6 on pg 306 (top of page)
- Q#4,5,7, 8ab on pg 306 (blue area)